




WE SERVE HUMANITY WITH EXCELLENCE

IBRAHIM MEDS



19- CALIPERS OTOLOGY

Calipers



IM-19-150-06
Castroviejo
8cm 3 1/4 "
Caliper measuring range
15 mm



IM-19-150-07
Castroviejo
8cm 4 1/4 "
Caliper measuring range
20 mm



IM-19-150-09
Castroviejo
8cm 5 1/4 "
Caliper measuring range
20 mm



IM-19-155-09
Castroviejo-Epker
9cm 3 1/2 "
Caliper measuring range
40 mm


Calipers




IM-19-155-18
Castroviejo-Epker
18cm 7 "
Caliper measuring range
40 mm




IM-19-157-18
Castroviejo-Epker
18cm 7 "
Caliper measuring range
40 mm




Otology




IM-19-160-01
Micro Ear Forceps
4,0 x 0,8 mm
smooth, straight
shaft length 80 mm




IM-19-160-06
Micro Ear Forceps
4,0 x 0,8 mm
serrated, straight
shaft length 80 mm




IM-19-160-07
Micro Ear Forceps
4,0 x 0,8 mm
bent upwards, serrated
shaft length 80 mm




IM-19-160-08
Micro Ear Forceps
4,0 x 0,8 mm
right, serrated
shaft length 80 mm




IM-19-160-09
Micro Ear Forceps
4,0 x 0,8 mm
left, serrated
shaft length 80 mm




IM-19-160-11
Mini-McGee
Micro Ear Forceps
3,5 x 0,6 mm
smooth, straight
shaft length 80 mm




IM-19-160-16
Mini-McGee
Micro Ear Forceps
3,5 x 0,6 mm
serrated, straight
shaft length 80 mm




IM-19-160-21
Micro Ear Forceps
4,0 x 0,6 mm
smooth, straight
shaft length 80 mm




IM-19-160-26
Micro Ear Forceps
4,0 x 0,6 mm
serrated, straight
shaft length 80 mm



IM-19-160-27
Micro Ear Forceps
3,0 x 0,6 mm
bent upwards, serrated
shaft length 80 mm




IM-19-160-28
Micro Ear Forceps
3,0 x 0,6 mm
right, serrated
shaft length 80 mm




IM-19-160-29
Micro Ear Forceps
3,0 x 0,6 mm
left, serrated
shaft length 80 mm




Otology




IM-19-165-01
Micro Ear Forceps
4,0 mm cup 1,0 x 0,9 mm
straight shaft length 80 mm




IM-19-165-02
Micro Ear Forceps
4,0 mm cup 1,0 x 0,9 mm
bent upwards
shaft length 80 mm




IM-19-165-03
Micro Ear Forceps
4,0 mm
cup 1,0 x 0,9 mm
right shaft length 80 mm




IM-19-165-04
Micro Ear Forceps
4,0 mm
cup 1,0 x 0,9 mm
left shaft length 80 mm




IM-19-165-11
Micro Ear Forceps
3,5 mm
cup 0,6 x 0,5 mm
straight shaft length 80 mm



IM-19-165-12
Micro Ear Forceps
3,5 mm
cup 0,6 x 0,5 mm
bent upwards
shaft length 80 mm




IM-19-165-13
Micro Ear Forceps
3,5 mm
cup 0,6 x 0,5 mm
right shaft length 80 mm




IM-19-165-14
Micro Ear Forceps
3,5 mm
cup 0,6 x 0,5 mm
left shaft length 80 mm




Otology




IM-19-170-01
Bellucci
Micro Ear Scissor
svery fine
4,0 x 0,8 mm
straight
shaft length 80 mm



IM-19-170-02
Bellucci
Micro Ear Scissor
svery fine
4,0 x 0,8 mm
bent upwards
shaft length 80 mm




IM-19-170-03
Bellucci
Micro Ear Scissor
svery fine
4,0 x 0,8 mm
right shaft
length 80 mm




IM-19-170-04
Bellucci
Micro Ear Scissor
svery fine
4,0 x 0,8 mm
left shaft
length 80 mm







IM-19-175-51
Wullstein
Micro Ear Scissors
4,0 x 1,0 mm
straight shaft 0,8 mm Ø
shaft length 80 mm



IM-19-175-52
Wullstein
Micro Ear Scissors
4,0 x 1,0 mm
bent upwards
shaft 0,8 mm Ø
shaft length 80 mm




IM-19-175-53
Wullstein
Micro Ear Scissors
4,0 x 1,0 mm
right shaft 0,8 mm Ø
shaft length 80 mm




IM-19-175-54
Wullstein
Micro Ear Scissors
4,0 x 1,0 mm
lefts haft 0,8 mm Ø
shaft length 80 mm





IM-19-180-11
McGee
Wire Bending Forceps
0,8 x 3,5 mm
bent downwards
Wire Bending Forceps
shaft length 80 mm



IM-19-180-12
McGee
Wire Bending Forceps
0,8 x 6,0 mm
bent downwards
Wire Bending Forceps
shaft length 80 mm



IM-19-185-01
House-Dieter
Malleus Nipper
up-cutting
Malleus Nipper
shaft length 80 mm



IM-19-185-02
House-Dieter
Malleus Nipper
down-cutting
Malleus Nipper
shaft length 80 mm



IM-19-185-06
House-Dieter
Malleus Nipper
right-cutting
Malleus Nipper
shaft length 80 mm



IM-19-185-07
House-Dieter
Malleus Nipper
left-cutting
Malleus Nipper
shaft length 80 mm





IM-19-185-01
House-Dieter
Malleus Nipper
up-cutting
Malleus Nipper
shaft length 80 mm



IM-19-185-02
House-Dieter
Malleus Nipper
down-cutting
Malleus Nipper
shaft length 80 mm



IM-19-185-06
House-Dieter
Malleus Nipper
right-cutting
Malleus Nipper
shaft length 80 mm



IM-19-185-07
House-Dieter
Malleus Nipper
left-cutting
Malleus Nipper
shaft length 80 mm



Otology



IM-19-190-02
Schuknecht
Rotary Knife
Rotary Knife



IM-19-190-06
Schuknecht
14.5cm 5 3/4"
Sickle Knife
Sickle Knife



IM-19-190-11
Wullstein
15.5cm 6"
Sickle Knife
Sickle Knife



IM-19-190-16
Tabb
Mucosa Knife
Mucosa Knife



IM-19-190-21
Plester
15cm 6"
Flap Knife
Flap Knife




ZI-22-022-01
Rosen
15cm 6"
Circular Cut
Knife
Fig. 13 mm Ø



ZI-22-022-02
Rosen
15cm 6"
Circular Cut
Knife
Fig. 2
2,8 mm Ø




ZI-22-022-03
Rosen
15cm 6"
Circular Cut
Knife
Fig. 3
2,6 mm Ø



ZI-22-022-04
Rosen
15cm 6"
Circular Cut
Knife
Fig. 4
2 mm Ø



IM-19-195-05
Rosen
15cm 6"
Circular Cut
Knife
Fig. 5
1,6 mm Ø



IM-19-195-06
Rosen
15cm 6"
Circular Cut
Knife
Fig. 6
1,4 mm Ø



IM-19-195-07
Rosen
15cm 6"
Circular Cut
Knife
Fig. 7
1,0 mm Ø

Otology



IM-19-200-01
Barbara
16cm 6 1/4"
Micro Ear Needle



IM-19-200-02
Barbara
16cm 6 1/4"
Micro Ear Needle



IM-19-200-03
Barbara
16cm 6 1/4"
Micro Ear Needle



IM-19-200-04
Rosen
16cm 6 1/4"
Micro Ear Needle



IM-19-200-06
Barbara
16cm 6 1/4"
Micro Ear Hook
0,3 mm



IM-19-200-07
Barbara
16cm 6 1/4"
Micro Ear Hook
0,6 mm



IM-19-200-08
Barbara
16cm 6 1/4"
Micro Ear Hook
1,0 mm



IM-19-200-11
Barbara
16cm 6 1/4"
Micro Ear Hook
0,3 mm



IM-19-200-12
Barbara
16cm 6 1/4"
Micro Ear Hook
0,6 mm



IM-19-200-13
Barbara
16cm 6 1/4"
Micro Ear Hook
1,0 mm



IM-19-200-14
Barbara
16cm 6 1/4"
Micro Ear Hook
1,2 mm



IM-19-200-15
Barbara
16cm 6 1/4"
Micro Ear Hook
1,5 mm



IM-19-200-16
Barbara
16cm 6 1/4"
Micro Ear Hook
2,0 mm



IM-19-200-17
Barbara
16cm 6 1/4"
Micro Ear Hook
2,5 mm

Otology



IM-19-200-01
Barbara
16cm 6 1/4"
Micro Ear Needle



IM-19-200-02
Barbara
16cm 6 1/4"
Micro Ear Needle



IM-19-200-03
Barbara
16cm 6 1/4"
Micro Ear Needle



IM-19-200-04
Rosen
16cm 6 1/4"
Micro Ear Needle



IM-19-200-06
Barbara
16cm 6 1/4"
Micro Ear Hook
0,3 mm



IM-19-200-07
Barbara
16cm 6 1/4"
Micro Ear Hook
0,6 mm



IM-19-200-08
Barbara
16cm 6 1/4"
Micro Ear Hook
1,0 mm



IM-19-200-11
Barbara
16cm 6 1/4"
Micro Ear Hook
0,3 mm



IM-19-200-12
Barbara
16cm 6 1/4"
Micro Ear Hook
0,6 mm



IM-19-200-13
Barbara
16cm 6 1/4"
Micro Ear Hook
1,0 mm



IM-19-200-14
Barbara
16cm 6 1/4"
Micro Ear Hook
1,2 mm




IM-19-200-15
Barbara
16cm 6 1/4"
Micro Ear Hook
1,5 mm





IM-19-200-16
Barbara
16cm 6 1/4"
Micro Ear Hook
2,0 mm



IM-19-200-17
Barbara
16cm 6 1/4"
Micro Ear Hook
2,5 mm







IM-19-210-31
Barbara
16cm 6 1/4"
Micro Ear Hook
0,5 mm

IM-19-210-32
Barbara
16cm 6 1/4"
Micro Ear Hook
1,0 mm



IM-19-215-01
15.5cm 6"
Micro Ear Hook
0,3 mm



IM-19-215-02
15.5cm 6"
Micro Ear Hook
0,6 mm



IM-19-215-03
15.5cm 6"
Micro Ear Hook
1,0 mm



IM-19-215-06
15.5cm 6"
Micro Ear Hook
0,3 mm



IM-19-215-07
15.5cm 6"
Micro Ear Hook
0,6 mm



IM-19-215-08
15.5cm 6"
Micro Ear Hook
1,0 mm



IM-19-220-01
Barbara
15.5cm 6"
Micro Ear Needle

IM-19-220-06
Barbara
15.5cm 6"
Micro Ear Needle



IM-19-225-06
Micro Ear Hook
0,3 mm 25°



IM-19-225-07
Micro Ear Hook
0,6 mm 25°



IM-19-225-08
Shambaugh
Micro Ear Hook
1,0 mm 25°



IM-19-225-11
Micro Ear Hook
0,3 mm 45°



IM-19-225-12
Micro Ear Hook
0,6 mm 45°



IM-19-225-13
Shambaugh
Micro Ear Hook
1,0 mm 45°

Otology



IM-19-230-20
Shea
Exploring Hook
2,0 mm 45°
blunt Exploring Hook



IM-19-230-30
Austin
Attic Dissector
3,0 mm 90°
blunt Attic Dissector



IM-19-235-01
Fisch
15.5cm 6"
Micro Dissector
angled to right



IM-19-235-02
Fisch
15.5cm 6"
Micro Dissector
angled to left



IM-19-235-05
Fisch
15.5cm 6"
Micro Dissector
angled to right



IM-19-235-06
Fisch
15.5cm 6"
Micro Dissector
angled to left



IM-19-240-23
Rosen
15.5cm 6"
Micro Ear Periosteal Elevator

Otology



IM-19-245-01
House
15.5cm 6"
Micro
Measuring Rod
3,5 mm



IM-19-245-02
House
15.5cm 6"
Micro
Measuring Rod
4,0 mm




IM-19-245-03
House
15.5cm 6"
Micro
Measuring Rod
4,5 mm



IM-19-245-04
House
15.5cm 6"
Micro
Measuring Rod
5,0 mm

Otology




IM-19-250-01
House
14.5cm 5 3/4"
Micro Curette
1,0 / 1,2 mm



IM-19-250-02
House
14.5cm 5 3/4"
Micro Curette
1,5 / 1,8 mm



IM-19-250-03
House
14.5cm 5 3/4"
Micro Curette
1,0 / 1,2 mm



IM-19-250-04
House
14.5cm 5 3/4"
Micro Curette
1,5 / 1,8 mm

Otology



IM-19-255-03
Toynbee
3.6cm 1 1/2"
Ear Specula,
Setfor adults
Fig. 2 – 4
Ear Specula, set



IM-19-255-11
Toynbee
3.6cm 1 1/2"
Ear Speculum
Fig. 14,0 mm Ø



IM-19-255-12
Toynbee
3.6cm 1 1/2"
Ear Speculum
Fig. 25,0 mm Ø



IM-19-255-13
Toynbee
3.6cm 1 1/2"
Ear Speculum
Fig. 36,0 mm Ø



IM-19-255-14
Toynbee
3.6cm 1 1/2"
Ear Speculum
Fig. 47,5 mm Ø



IM-19-255-04
Toynbee
3.6cm 1 1/2"
Ear Specula,
Setfor adults
Fig. 1 – 4
Ear Specula, set



IM-19-260-03
Hartmann
3.6cm 1 1/2"
Ear Specula,
Setfor children
Fig. 0 – 2
Ear Specula, set



IM-19-260-10
Hartmann
3.6cm 1 1/2"
Ear Speculum
Fig. 03,0 mm Ø



IM-19-260-11
Hartmann
3.6cm 1 1/2"
Ear Speculum
Fig. 14,0 mm Ø



IM-19-260-12
Hartmann
3.6cm 1 1/2"
Ear Speculum
Fig. 25,0 mm Ø



IM-19-260-13
Hartmann
3.6cm 1 1/2"
Ear Speculum
Fig. 36,0 mm Ø



IM-19-260-14
Hartmann
3.6cm 1 1/2"
Ear Speculum
Fig. 47,5 mm Ø



IM-19-260-03
Hartmann
3.61 1/2"
Ear Specula,
Setfor adults
Fig. 2 – 4
Ear Specula, set



IM-19-260-04
Hartmann
3.6cm 1 1/2"
Ear Specula,
Setfor adults
Fig. 1 – 4
Ear Specula, set



IM-19-265-00
Zöllner
3.8cm 1 1/2"
Ear Specula,
Setconsisting
of 6 pieces
Fig. 1 - 6
Ear Specula,
setoval, black



M-19-265-01
Zöllner
3.8cm 1 1/2"
Ear Speculum
Fig. 1
5,5 x 6,5 mm
oval, black



IM-19-265-02
Zöllner
3.8cm 1 1/2"
Ear Speculum
Fig. 2
6,5 x 7,5 mm
oval, black



IM-19-265-03
Zöllner
3.8cm 1 1/2"
Ear Speculum
Fig. 3
7,5 x 8,5 mm
oval, black



IM-19-265-04
Zöllner
3.8cm 1 1/2"
Ear Speculum
Fig. 48,5 x 9,5 mm
oval, black



IM-19-265-05
Zöllner
3.8cm 1 1/2"
Ear Speculum
Fig. 59,5 x 10,5 mm
oval, black



IM-19-265-06
Zöllner
3.8cm 1 1/2"
Ear Speculum
Fig. 6
10,5 x 11,5 mm
oval, black



1

IM-19-270-01
 Farrior
 3.6cm 1 1/2"
 Ear Speculum
 Fig. 15,5 mm
 black

2

IM-19-270-02
 Farrior
 3.6cm 1 1/2"
 Ear Speculum
 Fig. 2
 6,5 mm black

3

IM-19-270-03
 Farrior
 3.6cm 1 1/2"
 Ear Speculum
 Fig. 3
 7,5 mm black

4

IM-19-270-04
 Farrior
 3.6cm 1 1/2"
 Ear Speculum
 Fig. 4
 8,5 mm black

5

IM-19-270-05
 Farrior
 3.6cm 1 1/2"
 Ear Speculum
 Fig. 5
 9,5 mm black

6

IM-19-270-06
 Farrior
 3.6cm 1 1/2"
 Ear Speculum
 Fig. 6
 10,5 mm black

7

IM-19-270-07
 Farrior
 3.6cm 1 1/2"
 Ear Speculum
 Fig. 7
 11,5 mm black



IM-19-300-07
French
Tuning Fork
C2 - 16

IM-19-300-08
French
Tuning Fork
C1 - 32

IM-19-300-09
Hartmann
Tuning Fork
C 64

IM-19-300-10
Hartmann
Tuning Fork
c 128

IM-19-300-11
Hartmann
Tuning Fork
c1 - 256



IM-19-305-10
Martin
Tuning Fork
c 128



IM-19-305-11
Martin
Tuning Fork
c1 256



IM-19-305-12
Martin-Alloy
Tuning Fork
c2 512



IM-19-305-13
Martin-Alloy
Tuning Fork
c3 1024



IM-19-305-14
Martin-Alloy
Tuning Fork
c4 2048

Otology



IM-19-310-10
Hartmann
Tuning Fork
c 128



IM-19-310-11
Hartmann
Tuning Fork
c1 256



IM-19-310-12
Hartmann
Tuning Fork
c2 512



IM-19-315-13
Hartmann
Tuning Fork
c3 1024



IM-19-315-14
Hartmann
Tuning Fork
c4 2048



IM-19-315-15
Hartmann
Tuning Fork
c5 4096

Otology



IM-19-315-13
Hartmann
Tuning Fork
c3 1024



IM-19-315-14
Hartmann
Tuning Fork
c4 2048



IM-19-315-15
Hartmann
Tuning Fork
c5 4096



IM-15-320-10
Lucae
Tuning Fork
c 128



IM-19-320-11
Lucae
Tuning Fork
c 256



IM-19-320-12
Lucae
Tuning Fork
c 512



IM-19-320-13
Lucae
Tuning Fork
c 1024

Otology



IM-19-325-14
Lucae
Tuning Fork
c4 2048



IM-19-325-15
Lucae
Tuning Fork
c5 4096



IM-19-325-10
Lucae
Tuning Fork
c 128
movable up to h

Otology



IM-19-325-14
Lucae
Tuning Fork
c4 2048



IM-19-325-15
Lucae
Tuning Fork
c5 4096



IM-19-325-10
Lucae
Tuning Fork
c 128
movable up to h

Otology



IM-19-325-14
Lucae
Tuning Fork
c4 2048



IM-19-325-15
Lucae
Tuning Fork
c5 4096



IM-19-325-10
Lucae
Tuning Fork
c 128
movable up to h

Otology



IM-19-325-14
Lucae
Tuning Fork
c4 2048



IM-19-325-15
Lucae
Tuning Fork
c5 4096



IM-19-325-10
Lucae
Tuning Fork
c 128
movable up to h

Otology




IM-19-330-01
Hartmann
15cm 6"
Ear Catheter
Fig. 1



IM-19-330-02
Hartmann
15cm 6"
Ear Catheter
Fig. 2



IM-19-330-03
Hartmann
15cm 6"
Ear Catheter
Fig. 3





IM-19-335-16
Krause
18.5cm 7 1/4"
Polypus Snare
open tip



IM-19-336-16
Krause-Voss
18.5cm 7 1/4"
Polypus Snare
closed tip



Otology



IM-19-340-30
Snare Wire
0,3 mm \varnothing
coil of 10 m
Snare Wire



IM-19-342-40
Snare Wire
0,4 mm \varnothing
coil of 10 m
Snare Wire



IM-19-345-11
Troeltsch (Wilde)
10.5cm 4 1/8"
Ear Forceps



IM-19-345-12
Troeltsch (Wilde)
11cm 4 3/8"
Ear Forceps



IM-19-345-13
Troeltsch (Wilde)
12.5cm 5"
Ear Forceps



Otology



IM-19-350-14
Grünwald
14.5cm 5 3/4"
Ear Forceps



IM-19-355-16
Jansen
16cm 6 1/4"
Ear Forceps



IM-19-360-20
Jansen
21.5cm 8 1/2"
Ear Forceps



IM-19-365-12
Littauer
10cm 4"
Ear Polypus Forceps

Otology



IM-19-350-14
Grünwald
14.5cm 5 3/4"
Ear Forceps



IM-19-355-16
Jansen
16cm 6 1/4"
Ear Forceps



IM-19-360-20
Jansen
21.5cm 8 1/2"
Ear Forceps



IM-19-365-12
Littauer
10cm 4"
Ear Polypus Forceps




IM-19-370-12
Hartmann
12cm 4 3/4"
Ear Polypus Forceps
fine model

IM-19-370-14
Hartmann
14cm 5 1/2"
Ear Polypus Forceps
regular model



IM-19-375-11
Buck
11.5cm 4 1/2"
Foreign Body Lever
Foreign Body Forceps



IM-19-380-10
Quire
10.5cm 4 1/8"
Foreign Body Lever
Foreign Body Lever

Otology



IM-19-285-08
Hartmann-Wullstein
Ear Polypus Forceps
shaft length 85 mm



IM-19-286-08
Hartmann-Wullstein
Ear Polypus Forceps
shaft length
85 mm 1,0 mm Ø



IM-19-287-08
Hartmann
Ear Polypus Forceps
shaft length 85 mm



IM-19-288-14
Hartmann
Ear Polypus Forceps
shaft length
140 mm



IM-19-389-08
Hartmann
Ear Polypus Forceps
shaft length
85 mm 1:2 teeth



IM-19-390-81
Hartmann-Herzfeld
Ear Polypus Forceps
shaft length
85 mm 2,0 mm Ø



IM-19-391-82
Hartmann-Herzfeld
Ear Polypus Forceps
shaft length
85 mm 2,5 mm Ø



IM-19-392-83
Hartmann-Herzfeld
Ear Polypus Forceps
shaft length
85 mm 3,0 mm Ø



IM-19-393-08
Strümpel
Ear Polypus Forceps




IM-19-394-08
Hartmann-Hoffman
Ear Polypus Forceps
2,0 mm Ø
through-cutting





IM-19-400-17
Citelli
Mastoid Rongeur
3 mm cut



IM-19-405-19
Beyer
Mastoid Rongeur
2 mm cut

IM-19-405-14
Beyer
Mastoid Rongeur
1 mm cut



IM-19-410-13
Langen beck
13cm 5"
Ear Curetteblunt



IM-19-415-01
Shapleigh
15cm 6"
Ear Curette



IM-19-415-02
Shapleigh
15cm 6"
Ear Curette



IM-19-420-02
Buck
17xm 6 3/4"
Ear Curette
Fig. 00
blunt



IM-19-420-10
Buck
17cm 6 3/4"
Ear Curette
Fig. 0
blunt




IM-19-420-11
Buck
17cm 6 3/4"
Ear Curette
Fig. 1
blunt




IM-19-420-12
Buck
17cm 6 3/4"
Ear Curette
Fig. 2blunt


Otology


IM-19-325-01
Buck
17cm 6 3/4"
Ear Curette
Fig. 000blunt




IM-19-325-02
Buck
17cm 6 3/4"
Ear Curette
Fig. 00blunt




IM-19-325-10
Buck
17cm 6 3/4"
Ear Curette
Fig. 0blunt




IM-19-325-11
Buck
17cm 6 3/4"
Ear Curette
Fig. 1blunt




IM-19-425-12
Buck
17cm 6 3/4"
Ear Curette
Fig. 2blunt




IM-19-425-13
Buck
17cm 6 3/4"
Ear Curette
Fig. 3blunt




IM-19-425-14
Buck
17cm 6 3/4"
Ear Curette
Fig. 4blunt



IM-19-427-01
Buck
17cm 6 3/4"
Ear Curette
Fig. 000sharp




IM-19-427-02
Buck
17cm 6 3/4"
Ear Curette
Fig. 00sharp




IM-19-427-10
Buck
17cm 6 3/4"
Ear Curette
Fig. 0sharp



IM-19-427-11
Buc
k17cm 6 3/4"
Ear Curette
Fig. 1sharp



IM-19-427-12
Buck
17cm 6 3/4"
Ear Curette
Fig. 2sharp



IM-19-427-13
Buck
17cm 6 3/4"
Ear Curette
Fig. 3sharp



M-19-427-14
Buck
17cm 6 3/4"
Ear Curette
Fig. 4sharp

Otology



IM-19-430-15
Jansen
17cm 6 3/4"
Ear Scoopsharp



IM-19-435-12
Gross
13cm 5"
Ear Scoop/Hookblunt
Ear Scoop/Hook



IM-19-440-01
Weber-Loch
16.5cm 6 1/2"
Foreign Body
Instrument
Ear Loop



IM-19-445-01
Billeau
16cm 6 1/4"
Ear Loop
Ear Loop




IM-19-445-02
Billeau
16cm 6 1/4"
Ear Loop
Ear Loop



IM-19-445-03
Billeau
16cm 6 1/4"
Ear Loop
Ear Loop



Otology



IM-19-450-16
Moldenhauer
15.5cm 6"
Foreign Body Lever
Foreign Body Lever

IM-19-450-16
Zaufal
15cm 6"
Foreign Body Lever
Foreign Body Lever

IM-19-450-16
Hartmann
16cm 6 1/4"
Ear Hooklet
sharp Ear Hooklet




IM-19-455-01
Lucae
14cm 5 1/2"
Ear Hooklet small



IM-19-455-02
Lucae
14cm 5 1/2"
Ear Hooklet medium



IM-19-455-03
Lucae
14cm 5 1/2"
Ear Hooklet large



IM-19-460-01
Wagener
14cm 5 1/2"
Ear Hook let
Fig. 1

IM-19-460-02
Wagener
14cm 5 1/2"
Ear Hook let
Fig. 2

IM-19-460-03
Wagener
14cm 5 1/2"
Ear Hook let
Fig. 3


IM-19-460-04
Wagener
14cm 5 1/2"
Ear Hook let
Fig. 4

IM-19-460-05
Wagener
14cm 5 1/2"
Ear Hook let
Fig. 5

IM-19-460-06
Wagener
14cm 5 1/2"
Ear Hook let
Fig. 6 sharp

IM-19-460-07
Wagener
14cm 5 1/2"
Ear Hook let
Fig. 7

Otology




IM-19-465-01
Day
17cm 6 3/4"
Ear Hook let
long
IM-19-465-02
Day
17cm 6 3/4"
Ear Hook let short



IM-19-470-16
Poltzer
14.5cm 5 3/4"
Furuncle Knife
Furuncle Knife



IM-19-475-16
Poltzer
17.5cm 7"
Tympanum Lancet



IM-19-480-16
Poltzer
15.5cm 6"
Tympanum Needle
IM-19-482-16
Poltzer
16cm 6 1/4"
Tympanum Needle
needle horizontally

Otology



IM-19-485-18
Lucae
17.5cm 7"
Tympanum Needle



IM-19-487-18
Lucae
17.5cm 7"
Tympanum Needle

Otology



IM-19-490-20
Lempert
20.5cm 8"
Micro Ear Periosteal Elevator
6 mm




IM-19-490-17
Lempert
18cm 7"
Micro Ear Periosteal Elevator
6 mm

Otology



o	o	o	o	o
IM-19-495-10	IM-19-495-20	IM-19-495-30	IM-19-495-40	IM-19-495-50
Lempert	Lempert	Lempert	Lempert	Lempert
21.5cm 8 1/2"	21.5cm 8 1/2"	21.5cm 8 1/2"	21.5cm 8 1/2"	21.5cm 8 1/2"
Ear Scoop	Ear Scoop	Ear Scoop	Ear Scoop	Ear Scoop
sharp	sharp	sharp	sharp	sharp




IM-19-500-00
Ear Syringe
self fillin
gin plastic pouch



IM-19-505-00
Ear Syringe
complete with
protection disc
and 2 tips 50 ml

IM-19-505-04
Ear Syringe
complete with
protection disc
and 2 tips 100 ml

Otology



IM-19-510-01
Goldsmith
(Thudichum)
6.5cm 2 1/2"
Nasal Speculum
Fig. 1
blades
21,0 x 9,5 mm



IM-19-510-02
Goldsmith
(Thudichum)
6.5cm 2 1/2"
Nasal Speculum
Fig. 2
blades
21,0 x 10,5 mm



IM-19-510-03
Goldsmith
(Thudichum)
6.5cm 2 1/2"
Nasal Speculum
Fig. 3
blades
21,0 x 11,5 mm



IM-19-515-01
Hartmann-Halle
15.5cm 6"
Nasal Speculum
Fig. 1
length of blade
27 mm



IM-19-515-02
Hartmann-Halle
15.5cm 6"
Nasal Speculum
Fig. 2
length of blade
31 mm



IM-19-515-03
Hartmann-Halle
15.5cm 6"
Nasal Speculum
Fig. 3
length of blade
35 mm



Otology



IM-19-520-13
Martin
13cm 5"
Nasal Speculum
length of blade 17 mm




IM-19-522-13
Martin
12.5cm 5"
Nasal Speculum
length of blade 17 mm



IM-19-524-13
Tieck-Halle
13cm 5"
Nasal Speculum
length of blade 16 mm



IM-19-525-01
Killian
13cm 5"
Nasal Speculum
Fig. 1
length of blade
35 mm



IM-19-526-02
Killian
13cm 5"
Nasal Speculum
Fig. 2
length of blade
50 mm



IM-19-527-03
Killian
13cm 5"
Nasal Speculum
Fig. 3
length of blade
75 mm




IM-19-528-04
Killian
13cm 5"
Nasal Speculum
Fig. 4
length of blade
85 mm



IM-19-530-14
Modell Wien
13.5cm 5 1/4"
Nasal Speculum
length of blade 30 mm

Otology



IM-19-535-35
Cottle
13.5cm 5 1/4"
Nasal Speculum
Fig. 1
length of blade
35 mm

IM-19-535-50
Cottle
13.5cm 5 1/4"
Nasal Speculum
Fig. 2
length of blade
50 mm

IM-19-535-75
Cottle
13.5cm 5 1/4"
Nasal Speculum
Fig. 3
length of blade
75 mm

IM-19-535-90
Cottle
14.5cm 5 3/4"
Nasal Speculum
Fig. 4
length of blade
85 mm





IM-19-540-45
Cottle
13.5cm 5 1/4"
length of blade
45 mm



IM-19-540-75
Cottle
14cm 5 1/2"
length of blade
75 mm

Otology




IM-19-545-26
Krause
26cm 10 1/4"
Polypus Snare
shaft 3 mm Ø

IM-19-550-26
Krause-Voss
26cm 10 1/4"
Polypus Snare
shaft 3 mm Ø



IM-19-552-12
10.5cm 4 1/8"
Snare Wire
0,3 mm Ø
bundle with
12 pieces
Snare Wire



IBRAHIM MEDS

I N S T R U M E N T S

 77 A Model Town Ugoki Sialkot

 +92 52 3512235

 +92 52 4007564

 +92 300 7163508 , +92 307 6427547

 info@ibrahimmeds.com

 amjad@ibrahimmeds.com

 www.ibrahimmeds.com